

CAMARA DE PRUEBA DE ALTA-BAJA TEMPERATURA Y BAJA PRESION DE AIRE SERIE TET4 Y TET7



No.	Product name	Type	Dimensions of work chamber (mm)	Power (about Kw)	Temperature range (°C)	Other parameters	Implementation / Satisfaction Standard
1	High-low temperature and low air pressure test chamber	TET401Y	450×450×580	6.5	-40~+130	Extreme pressure 0.5Kpa	GB/T2423.25-92 GB/T2423.26-92 GJB150.2-86
2		TET402Y	500×600×650	7.5			
3		TET405Y	900×750×750	12	-40~+150		
4		TET410Y	1000×1000×1000	16	-70~+130	Equipped with high-efficiency vacuum pump with well-known brand	
5		TET701Y	450×450×580	8.5			
6		TET702Y	500×600×650	10	-70~+150		
7		TET705Y	900×750×750	14.5			
8		TET710Y	1000×1000×1000	19			