



AGITADOR DIGITAL PARA TAMICES DE 8", ELECTROMAGNETICO CON MOVIMIENTO EN 3 EJES,
110/120 VOLTS 60HZ CATALOGO OCT11060 MARCA ENDECOTTS MODELO OCTAGON 200



Octagon 200

- Range: 20 μ m to 125 mm
- Drive / sieving motion: electromagnetic 3D
- Amplitude: 0 - 3 mm, digital setting in 10 steps
- Sieve diameter: 100 / 200 mm, 3" / 8"